

Fastest

Simple
20 minute
Midi/Maxi
preps

Highest Yield

6x more
plasmid

Ultra-Pure

Vaccine grade
& transfection
ready

* Endotoxin levels are typically <0.025 EU/μg DNA

ZymoPURE™ II	Binding capacity	Elution volume	Volume of <i>E.coli</i> culture	Process time	Size	Cat No.	2018 Special Price (Ex. GST)
Midiprep Kit	300μg	100-200 μl	50 ml	18 min.	25 Preps	D4200	\$308
					50 Preps	D4201	\$554
Maxiprep Kit	1200μg	200-400 μl	150 ml	18 min.	10 Preps	D4202	\$249
					20 Preps	D4203	\$482
Gigaprep Kit	10mg	2-5 ml	2.5 L	40-50 min.	5 Preps	D4204	\$559

Request your
Midi or Maxi
Prep sample
today!

Buy any ZymoPure™ II Kit
and receive a
15% discount
on the EZ-Vac™
Vacuum Manifold

Cat#S7000 Ez-Vac™ Vacuum Manifold
Discount price \$347 (Ex-GST)



Valid until 31 December 2018. Please quote CS#2064
on your order to receive this pricing

ZymoPURE™ II Plasmid Kits, Plasmid Purification Reinvented